

HI6221 Advanced pH & ORP Meter

Dear Customer, thank you for choosing Hanna Instruments.

For more information about Hanna Instruments and our products, visit www.hannainst.com or e-mail us at sales@hannainst.com. For technical support, contact your local Hanna Instruments office or e-mail us at tech@hannainst.com.

Please scan the QR code or use the link below to download the user manual.




<https://manuals.hannainst.com/HI6221>



Package Contents

Each HI6221 is supplied with:

- HI1131B pH electrode
- HI7662-TW temperature probe
- pH 4.01 buffer solution (2 sachets)
- pH 7.01 buffer solution (4 sachets)
- pH 10.01 buffer solution (2 sachets)
- Electrode cleaning solution (2 sachets)
- Electrode storage solution (25 mL)
- 3.5M KCl electrolyte fill solution (30 mL)
- Capillary pipette (1 pc.)
- Electrode arm with base plate and screw (requires mounting)
- 24 VDC power adapter
- USB-C to USB-A cable
- Probe quality certificates
- Instrument quality certificate
- Quick reference guide




Key Features

- Application profiles allow for quick and direct measurement without updating the settings
- Log files can be transferred to a USB flash drive
- Ethernet and Wi-Fi[®] network data sharing options

Getting Started


i. Identify the base plate (1), place the meter face down, and proceed as illustrated here.



ii. Follow the user manual for electrode preparation.

iii. Connect the pH electrode (2), temperature probe (3), and power adapter (4) to the rear panel of the meter.

iv. Connect the power plug and press the black ON/OFF power button (5).



Extension port

Note: Remove the transparent film that protects the capacitive touchscreen prior to operating the meter.

*This device, when running firmware version v1.3.0, is not equipped with Wi-Fi functionality. The Wi-Fi icon (📶) will not be present.

General Description

The HI6221 is Hanna Instruments advanced pH and ORP meter with a large touchscreen display and streamlined design.

System Overview

Hardware Flexibility

Connectors for printer, keyboard, Ethernet

User Interface

7-inch color display, 800 × 480 resolution

Multi-touch screen supports video playback and data plotting

Measurement Flexibility

User profiles

Accessible logging control

Select Operating Language

At first power on the meter launches into a startup tutorial.

By default English is selected.

Use the language window to select the operating language.

Touchscreen

Icon	Name	Function
	Back	returns user to previous hierarchical menu level exit or escape function
	Home	access to measurement screen exit or escape function
	Menu	access Users, System Settings, Log Recall, Reports, Help

Basic Operations

General operating modes are setup, measurement, logging, data sharing.

- Tap  key to access:
 -  User settings
 -  System settings
 -  Log recall files and file management.
 - User can view a single sample or an interval log session.
 -  Reports for method specific applications
 -  Help for text and video support
- Tap  key to return to measurement.
- Tap  icon (Measurement Settings) to access sensor-related functions.

Log in & Create New Account

- Tap  followed by .
- Tap **Edit Users** to enter Account Management.
- Tap  to enable Account Creation and Logins. Tap  to return.
- Tap **Logout** to enter Users screen.
 - “Admin” account is automatically created.
- Tap the **plus** symbol avatar.
- Input user name and tap .
- Enter password and tap . Re-enter password to confirm.

Electrode Usage Guidelines

- Remove protective cap and rinse with water to eliminate salt build up.
- Check correct electrolyte level.
 - Refill when level falls below 1/2” from the fill hole opening. Keep the fill hole open (1) during calibration and measurement.
- To refresh the electrode, soak in **HI70300** Storage solution for at least 30 minutes.
- Keep the glass bulb (3) and junction (2) hydrated.
 - Store with a few drops of storage solution in the protective cap.
- Keep the junction (2) immersed during calibration/measurement.



KEEP DEVICE UP TO DATE

To introduce new features and/or performance improvements Hanna Instruments releases updated firmware versions.

To check for new releases, scan the QR code or go to: <https://software.hannainst.com>



Hanna Instruments is committed to developing and deploying digital solutions with a positive impact on the environment and climate.

All Hanna instruments conform to the CE European Directives and our production facilities are ISO 9001 certified. HI6221 is warranted for a period of two years (probes for six months) against defects in workmanship and materials when used for its intended purpose and maintained according to instructions. 

Please retain for future use.

QR6221 08/25